



PAEP Antibody, Biotin conjugated

Product Code	CSB-PA017381LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09466
Immunogen	Recombinant Human Glycodelin protein (19-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycodelin (GD) (Placental protein 14) (PP14) (Pregnancy-associated endometrial alpha-2 globulin) (PAEG) (PEG) (Progestagen-associated endometrial protein) (Progesterone-associated endometrial protein) (Zona-binding inhibitory factor-1) (ZIF-1), PAEP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PAEP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.