



PDE2A Antibody, HRP conjugated

Product Code	CSB-PA017687LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00408
Immunogen	Recombinant Human cGMP-dependent 3',5'-cyclic phosphodiesterase protein (15-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cGMP-dependent 3',5'-cyclic phosphodiesterase (EC 3.1.4.17) (Cyclic GMP-stimulated phosphodiesterase) (CGS-PDE) (cGSPDE), PDE2A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE2A