



PEPD Antibody, HRP conjugated

Product Code	CSB-PA017784LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12955
Immunogen	Recombinant Human Xaa-Pro dipeptidase protein (258-395AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Xaa-Pro dipeptidase (X-Pro dipeptidase) (EC 3.4.13.9) (Imidodipeptidase) (Peptidase D) (Proline dipeptidase) (Prolidase), PEPD, PRD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEPD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.