





PLA2G10 Antibody

Product Code	CSB-PA018085LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15496
Immunogen	Recombinant Human Group 10 secretory phospholipase A2 protein (43-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Group 10 secretory phospholipase A2 (EC 3.1.1.4) (Group X secretory phospholipase A2) (GX sPLA2) (sPLA2-X) (Phosphatidylcholine 2-acylhydrolase 10), PLA2G10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLA2G10

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA018085LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA018085LA01HU at dilution of 1:100

