





PLD3 Antibody

Product Code	CSB-PA018146LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IV08
Immunogen	Recombinant Human Phospholipase D3 protein (70-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase D3 (PLD 3) (EC 3.1.4.4) (Choline phosphatase 3) (HindIII K4L homolog) (Hu-K4) (Phosphatidylcholine-hydrolyzing phospholipase D3), PLD3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLD3

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA018146LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA018146LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)