



PLD3 Antibody, HRP conjugated

Product Code	CSB-PA018146LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IV08
Immunogen	Recombinant Human Phospholipase D3 protein (70-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase D3 (PLD 3) (EC 3.1.4.4) (Choline phosphatase 3) (HindIII K4L homolog) (Hu-K4) (Phosphatidylcholine-hydrolyzing phospholipase D3), PLD3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.