



PNLIPRP1 Antibody, HRP conjugated

Product Code	CSB-PA018263LB01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06857
Immunogen	Recombinant Dog Inactive pancreatic lipase-related protein 1 protein (18-467AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive pancreatic lipase-related protein 1 (PL-RP1), PNLIPRP1, PLRP1
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others
Target Names	PNLIPRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.