



PNLIPRP1 Antibody, Biotin conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA018263LD01DO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P06857 |
| Immunogen | Recombinant Dog Inactive pancreatic lipase-related protein 1 protein (18-467AA) |
| Raised In | Rabbit |
| Species Reactivity | Dog |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Inactive pancreatic lipase-related protein 1 (PL-RP1), PNLIPRP1, PLRP1 |
| Immunogen Species | Canis lupus familiaris (Dog) (Canis familiaris) |
| Research Area | Others |
| Target Names | PNLIPRP1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |