




PNMA2 Antibody

Product Code	CSB-PA018267LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL42
Immunogen	Recombinant Human Paraneoplastic antigen Ma2 protein (2-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeuroal antigen Ma2) (Paraneoplastic neuronal antigen MM2), PNMA2, KIAA0883 MA2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PNMA2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA018267LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.