



PNMA2 Antibody, Biotin conjugated

Product Code	CSB-PA018267LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL42
Immunogen	Recombinant Human Paraneoplastic antigen Ma2 protein (2-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeuroal antigen Ma2) (Paraneoplastic neuronal antigen MM2), PNMA2, KIAA0883 MA2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PNMA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.