



PRL Antibody, HRP conjugated

Product Code	CSB-PA018724YB01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12420
Immunogen	Recombinant Horse Prolactin protein (31-229AA)
Raised In	Rabbit
Species Reactivity	Horse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin (PRL), PRL
Immunogen Species	Equus caballus (Horse)
Research Area	Others
Target Names	PRL