



PRL Antibody, Biotin conjugated

Product Code	CSB-PA018724YD01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12420
Immunogen	Recombinant Horse Prolactin protein (31-229AA)
Raised In	Rabbit
Species Reactivity	Horse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin (PRL), PRL
Immunogen Species	Equus caballus (Horse)
Research Area	Others
Target Names	PRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.