



PROC Antibody, HRP conjugated

Product Code	CSB-PA018742LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04070
Immunogen	Recombinant Human Vitamin K-dependent protein C protein (288-409AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitamin K-dependent protein C (EC 3.4.21.69) (Anticoagulant protein C) (Autoprothrombin IIA) (Blood coagulation factor XIV) [Cleaved into: Vitamin K-dependent protein C light chain; Vitamin K-dependent protein C heavy chain; Activation peptide], PROC
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PROC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.