



PRODH Antibody, HRP conjugated

Product Code	CSB-PA018744LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43272
Immunogen	Recombinant Human Proline dehydrogenase 1, mitochondrial protein (71-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proline dehydrogenase 1, mitochondrial (EC 1.5.5.2) (Proline oxidase) (Proline oxidase 2) (p53-induced gene 6 protein), PRODH, PIG6 POX2 PRODH2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRODH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.