



PROP1 Antibody, Biotin conjugated

Product Code	CSB-PA018753LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75360
Immunogen	Recombinant Human Homeobox protein prophet of Pit-1 protein (6-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein prophet of Pit-1 (PROP-1) (Pituitary-specific homeodomain factor), PROP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PROP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.