



PSMA2 Antibody, HRP conjugated

Product Code	CSB-PA018866LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25787
Immunogen	Recombinant Human Proteasome subunit alpha type-2 protein (120-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-2 (EC 3.4.25.1) (Macropain subunit C3) (Multicatalytic endopeptidase complex subunit C3) (Proteasome component C3), PSMA2, HC3 PSC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.