



PSMA2 Antibody, Biotin conjugated

Product Code	CSB-PA018866LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25787
Immunogen	Recombinant Human Proteasome subunit alpha type-2 protein (120-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-2 (EC 3.4.25.1) (Macropain subunit C3) (Multicatalytic endopeptidase complex subunit C3) (Proteasome component C3), PSMA2, HC3 PSC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA2