



PSMA3 Antibody, FITC conjugated

Product Code	CSB-PA018867LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25788
Immunogen	Recombinant Human Proteasome subunit alpha type-3 protein (218-248AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-3 (EC 3.4.25.1) (Macropain subunit C8) (Multicatalytic endopeptidase complex subunit C8) (Proteasome component C8), PSMA3, HC8 PSC8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA3