



PSMA3 Antibody, Biotin conjugated

Product Code	CSB-PA018867LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25788
Immunogen	Recombinant Human Proteasome subunit alpha type-3 protein (218-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-3 (EC 3.4.25.1) (Macropain subunit C8) (Multicatalytic endopeptidase complex subunit C8) (Proteasome component C8), PSMA3, HC8 PSC8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.