





PSMA4 Antibody

Product Code	CSB-PA018869LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25789
Immunogen	Recombinant Human Proteasome subunit alpha type-4 protein (148-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-4 (EC 3.4.25.1) (Macropain subunit C9) (Multicatalytic endopeptidase complex subunit C9) (Proteasome component C9) (Proteasome subunit L), PSMA4, HC9 PSC9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMA4

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA018869LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA018869LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.