



PSMD9 Antibody, Biotin conjugated

Product Code	CSB-PA018914LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00233
Immunogen	Recombinant Human 26S proteasome non-ATPase regulatory subunit 9 protein (103-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 9 (26S proteasome regulatory subunit p27), PSMD9
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.