





JAG1 Antibody

Product Code	CSB-PA01949A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78504
Immunogen	Recombinant Human Protein jagged-1 protein (809-1046AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein jagged-1 (Jagged1) (hJ1) (CD antigen CD339), JAG1, JAGL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	JAG1


Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA01949A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA01949A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA01949A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)