




RND2 Antibody

Product Code	CSB-PA019812LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52198
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoN protein (104-163AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoN (Rho family GTPase 2) (Rho-related GTP-binding protein Rho7) (Rnd2), RND2, ARHN RHO7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RND2

Image



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse liver tissue, Mouse brain tissue

All lanes: RND2 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.