



RND2 Antibody, HRP conjugated

Product Code	CSB-PA019812LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52198
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoN protein (104-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoN (Rho family GTPase 2) (Rho-related GTP-binding protein Rho7) (Rnd2), RND2, ARHN RHO7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RND2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.