



RNF126 Antibody, Biotin conjugated

Product Code	CSB-PA019828LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV68
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF126 protein (39-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF126 (EC 2.3.2.27) (RING finger protein 126), RNF126
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF126
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.