




RNF40 Antibody

Product Code	CSB-PA019890LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75150
Immunogen	Recombinant Human E3 ubiquitin-protein ligase BRE1B protein (471-724AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase BRE1B (BRE1-B) (EC 2.3.2.27) (95 kDa retinoblastoma-associated protein) (RBP95) (RING finger protein 40) (RING-type E3 ubiquitin transferase BRE1B), RNF40, BRE1B KIAA0661
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF40

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA019890LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.