



RNF40 Antibody, Biotin conjugated

Product Code	CSB-PA019890LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75150
Immunogen	Recombinant Human E3 ubiquitin-protein ligase BRE1B protein (471-724AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase BRE1B (BRE1-B) (EC 2.3.2.27) (95 kDa retinoblastoma-associated protein) (RBP95) (RING finger protein 40) (RING-type E3 ubiquitin transferase BRE1B), RNF40, BRE1B KIAA0661
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.