



RPL4 Antibody, Biotin conjugated

Product Code	CSB-PA020289LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36578
Immunogen	Recombinant Human 60S ribosomal protein L4 protein (12-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L4 (60S ribosomal protein L1) (Large ribosomal subunit protein uL4), RPL4, RPL1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.