




RPS15 Antibody


Product Code	CSB-PA020372LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62841
Immunogen	Recombinant Human 40S ribosomal protein S15 protein (63-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S15 (RIG protein) (Small ribosomal subunit protein uS19), RPS15, RIG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS15

Image




Western Blot


Positive WB detected in: HeLa whole cell lysate
 All lanes: RPS15 antibody at 4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 18 kDa
 Observed band size: 18 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA020372LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA020372LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA020372LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.