



RPS2 Antibody, Biotin conjugated

Product Code	CSB-PA020393LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15880
Immunogen	Recombinant Human 40S ribosomal protein S2 protein (2-61AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S2 (40S ribosomal protein S4) (Protein LLRep3) (Small ribosomal subunit protein uS5), RPS2, RPS4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RPS2