



SDCBP Antibody, Biotin conjugated

Product Code	CSB-PA020893LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00560
Immunogen	Recombinant Human Syntenin-1 protein (21-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntenin-1 (Melanoma differentiation-associated protein 9) (MDA-9) (Pro-TGF-alpha cytoplasmic domain-interacting protein 18) (TACIP18) (Scaffold protein Pbp1) (Syndecan-binding protein 1), SDCBP, MDA9 SYCL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SDCBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.