



SERPING1 Antibody, HRP conjugated

Product Code	CSB-PA021086LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05155
Immunogen	Recombinant Human Plasma protease C1 inhibitor protein (23-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma protease C1 inhibitor (C1 Inh) (C1Inh) (C1 esterase inhibitor) (C1-inhibiting factor) (Serp1n G1), SERPING1, C1IN C1NH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPING1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.