



SERPING1 Antibody, HRP conjugated

Product Code	CSB-PA021086LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05155
Immunogen	Recombinant Human Plasma protease C1 inhibitor protein (23-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma protease C1 inhibitor (C1 Inh) (C1Inh) (C1 esterase inhibitor) (C1-inhibiting factor) (Serp G1), SERPING1, C1IN C1NH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPING1