



SLC10A3 Antibody, Biotin conjugated

Product Code	CSB-PA021375LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09131
Immunogen	Recombinant Human P3 protein (7-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P3 protein (Solute carrier family 10 member 3), SLC10A3, DXS253E P3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC10A3