



SLC16A3 Antibody, FITC conjugated

Product Code	CSB-PA021410LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15427
Immunogen	Recombinant Human Monocarboxylate transporter 4 protein (406-465AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 4 (MCT 4) (Solute carrier family 16 member 3), SLC16A3, MCT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A3