



SLC16A7 Antibody, HRP conjugated

Product Code	CSB-PA021414LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60669
Immunogen	Recombinant Human Monocarboxylate transporter 2 protein (197-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 2 (MCT 2) (Solute carrier family 16 member 7), SLC16A7, MCT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.