



SLC16A7 Antibody, Biotin conjugated

Product Code	CSB-PA021414LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60669
Immunogen	Recombinant Human Monocarboxylate transporter 2 protein (197-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 2 (MCT 2) (Solute carrier family 16 member 7), SLC16A7, MCT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC16A7