



SLC18A2 Antibody, Biotin conjugated

Product Code	CSB-PA021427LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05940
Immunogen	Recombinant Human Synaptic vesicular amine transporter protein (463-514AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicular amine transporter (Monoamine transporter) (Solute carrier family 18 member 2) (Vesicular amine transporter 2) (VAT2), SLC18A2, SVMT VMAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC18A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.