




SLU7 Antibody

Product Code	CSB-PA021783LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95391
Immunogen	Recombinant Human Pre-mRNA-splicing factor SLU7 protein (138-310AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-splicing factor SLU7 (hSlu7), SLU7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLU7

Image



Western Blot

Positive WB detected in: K562 whole cell lysate, Mouse brain tissue


All lanes: SLU7 antibody at 3μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69 kDa


Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA021783LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA021783LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA021783LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.