





SOX4 Antibody

Product Code	CSB-PA022431LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06945
Immunogen	Recombinant Human Transcription factor SOX-4 protein (1-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor SOX-4, SOX4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SOX4

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA022431LA01HU at dilution of 1:100



Western Blot

Positive WB detected in: HeLa whole cell lysate, 293 whole cell lysate, HepG2 whole cell lysate
All lanes: SOX4 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 48 kDa
Observed band size: 48 kDa



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.