



SPRR1A Antibody, FITC conjugated

Product Code	CSB-PA022610LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35321
Immunogen	Recombinant Human Cornifin-A protein (8-72AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cornifin-A (19 kDa pancornulin) (SPRK) (Small proline-rich protein IA) (SPR-IA), SPRR1A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SPRR1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.