



STAT5B Antibody, FITC conjugated

Product Code	CSB-PA022815LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51692
Immunogen	Recombinant Human Signal transducer and activator of transcription 5B protein (40-280AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal transducer and activator of transcription 5B antibody; STA5B_HUMAN antibody; STAT5 antibody; Stat5b antibody; Transcription factor STAT5B antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAT5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.