



TERT Antibody, FITC conjugated

Product Code	CSB-PA023391LC01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6A548
Immunogen	Recombinant Dog Telomerase reverse transcriptase protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Dog
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Telomerase reverse transcriptase (EC 2.7.7.49) (Telomerase catalytic subunit), TERT
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others
Target Names	TERT