



TERT Antibody, Biotin conjugated

Product Code	CSB-PA023391LD01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6A548
Immunogen	Recombinant Dog Telomerase reverse transcriptase protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Telomerase reverse transcriptase (EC 2.7.7.49) (Telomerase catalytic subunit), TERT
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Others
Target Names	TERT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.