





TFEB Antibody

Product Code	CSB-PA023429EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19484
Immunogen	Recombinant Human Transcription factor EB protein (1-476AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor EB (Class E basic helix-loop-helix protein 35) (bHLHe35), TFEB, BHLHE35
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TFEB


Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA023429EA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA023429EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in: Rat lung tissue

All lanes: TFEB antibody at 3 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53, 44 kDa

Observed band size: 70 kDa