



TFEB Antibody, FITC conjugated

Product Code	CSB-PA023429EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19484
Immunogen	Recombinant Human Transcription factor EB protein (1-476AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor EB (Class E basic helix-loop-helix protein 35) (bHLHe35), TFEB, BHLHE35
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TFEB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.