



TGDS Antibody, FITC conjugated

Product Code	CSB-PA023444LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95455
Immunogen	Recombinant Human dTDP-D-glucose 4,6-dehydratase protein (143-305AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	dTDP-D-glucose 4,6-dehydratase (EC 4.2.1.46), TGDS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TGDS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.