



TPM3 Antibody, Biotin conjugated

Product Code	CSB-PA024106LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06753
Immunogen	Recombinant Human Tropomyosin alpha-3 chain protein (18-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tropomyosin alpha-3 chain (Gamma-tropomyosin) (Tropomyosin-3) (Tropomyosin-5) (hTM5), TPM3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TPM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.