



TPST1 Antibody, Biotin conjugated

Product Code	CSB-PA024132LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60507
Immunogen	Recombinant Human Protein-tyrosine sulfotransferase 1 protein (283-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-tyrosine sulfotransferase 1 (EC 2.8.2.20) (Tyrosylprotein sulfotransferase 1) (TPST-1), TPST1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TPST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.