



TYRO3 Antibody, HRP conjugated

Product Code	CSB-PA025396LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06418
Immunogen	Recombinant Human Tyrosine-protein kinase receptor TYRO3 protein (728-873AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase receptor TYRO3 (EC 2.7.10.1) (Tyrosine-protein kinase BYK) (Tyrosine-protein kinase DTK) (Tyrosine-protein kinase RSE) (Tyrosine-protein kinase SKY) (Tyrosine-protein kinase TIF), TYRO3, BYK DTK RSE SKY TIF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TYRO3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.